

DNA Isolation from ES Cells

From 96-well plate

- 1) Grow ES cells in a 96-well plate to be over-confluent. You don't need LIF in your media.
- 2) Quickly invert the plate over to dump media.
- 3) Wash cells 1x with ~ 150 uL PBS, then dump PBS.
- 4) Add 50 ul of Bradley Lysis Buffer containing proteinase K.
- 5) Replace lid and seal the plate with parafilm. Put the plate into a humidified chamber
- 6) Incubate in the humidified chamber overnight at 60 °C.
- 7) Allow the plate to cool to RT.
- 8) Add 100 uL ice-cold Ethanol/NaCl mix to precipitate DNA and mix well. Then incubate the plate about 30 minutes at RT.
- 9) Spin in 96-well plate holder centrifuge, 3000rpm 20 minutes.
- 10) Invert the plate to decant liquid. Blot on paper towels.
- 11) Add 150uL cold 70% EtOH and spin 10 minutes at 3000 rpm to rinse the pellet. Decant supernatant.
- 12) Repeat washing step (#11) and air-dry DNA pellet until there is no detectable Ethanol smell (approx 10-15 minutes)
- 13) Add 30 uL of warm TE pH 8.0
- 14) Cover Plate with wax cover and incubate at 56 °C ~10 mins

Bradley Lysis Buffer

<u>Stock Solutions</u>	<u>Final Conc</u>	<u>for 50mL</u>
1M Tris-HCl (pH 7.5)	10 mM	500 uL
0.5M EDTA	10 mM	1 mL
10% SDS	0.5 %	2.5 mL
5M NaCl	10 mM	100 uL
H2O		45.9 mL

ProteinaseK (20mg/mL)

:Add just before use as 1mg/ml final concentration

EtOH/NaCl Mix

100% EtOH	394mL
<u>5M NaCl</u>	<u>6 mL</u>
Total volume	400 mL